IMV China


new product  
introduction

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**Marshall Ma**  
Marketing Manager IMV China

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
2- Backup boar fertility prediction software—Ecotext

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3- Intermuscular fat measurement software—Biosoft

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

# WHAT parameters of semen analysis? 精液分析有哪些参数?

I want only  
“round strawberry flavored  
candy”!  
我只想要“圆的草莓味的糖果”!

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I want only  
“round strawberry  
flavored candy”!  
我只想要“圆的草莓味的糖果”!

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
I want only  
“round strawberry flavored  
candy”!  
我只想要“圆的草莓味的糖果”!

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



EK5 Viability & Acrosome integrity 活率和顶体完整性

EK6 Bacterial concentration 细菌浓度

# IMV strenghtening position with easyCyte

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




**Bovine牛:**  $R^2= 0.84$

*with high green fluorescence representative of immature cells). A significant relationship ( $R^2 = 0.84, P < 0.05$ ) was observed between real and predicted fertility. Once the accuracy of fertility prediction has been confirmed, the model developed in the present study could be used by artificial insemination centers for bull selection or for elimination of poor fertility ejaculates.*

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
**Equine马:**  $R^2= 0.942$

In conclusion, the results of the present experiment has enabled the determination of a new protocol for the evaluation of stallion semen, combining microscopical observation, computer-assisted motility analysis and flow cytometry, and providing a high level of fertility prediction. Similar investigations are in progress on frozen semen.

(Barrier B. et al. Theriogenology. 2016)

**Fig. 3.** Relation between the observed field fertility value (“observed fertility”) and the predicted fertility value based on the linear regression model (“calculated fertility”). For each stallion ( $n = 33$ ), the observed fertility is presented as a function of the calculated fertility, as well as the minimal and maximal values of the 90% confidence interval (“CI 90% mini” and “CI 90% maxi”). The equation of the regression curve between the observed fertility and the calculated fertility is presented with its  $r^2$ . The adjusted  $r^2$  of the linear regression model that was used to calculate the fertility was 0.942.

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
The newly developed Fertiboar technology is based on the analysis of boar testicular ultrasound pictures which reflect characteristics of the sperm producing tissue (Theriogenology 158 (2020): pp 58-65). A machine learning approach (Convolutional Neural Network) was used to identify regions within the testicular pictures that are key to semen production. Those regions were analysed to calculate parameters describing novel attributes like the homogeneity and echogenicity of the tissue.

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### MARBLING STANDARDS<sup>2</sup>


Courtesy of the National Pork Board, Des Moines, IA.

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Source: french alibaba 2017

BioSoft Toolbox® for Swine 4.0.1.1 © 2007-2020 by Biotronics, Inc.

File Tools Options Help

(N) Next Animal

Image Capture Tools  
(S) Stream (F) Freeze  
(C) Capture Image

Scanner Setup  
Depth: 14cm  
Probe Size: 13cm

Image Interpretation Tools  
 (I) IMF  
# of ROI: 2  
ROI Size: 66  
 w/ Depth Assist  
 (R) Area  
 (D) Fat Depth  
(A) Assist With Depths

Sort By Scan Time

Animal List  
16-4  
23-12  
24-9  
24-10  
27-1  
28-1  
28-6  
28-9  
28-10  
29-1

Animal Information  
ID: 28-10  
IMF: 0.3  
Area: 6.83  
Fat Depth: 0.49  
Loin Depth:  
Weight: 232  
Sex: 2  
Breed: Berkshire  
Farrow Date:  
Sire:  
Comments:

14-Oct-20 10:57 AM  
BioSoft Toolbox®: 4.0.1.0  
Biotronics DDI for Exago, a licensed product  
Exago Software Version: 01-12

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# BioSoft Toolbox<sup>®</sup> 4.0

**Biotronics** Inc.  
Ultrasound Solutions

## **B** Depth Assist Settings

### Measurement Locations

#### **Last Rib - Depth measurement taken anterior of the last rib**

- Between second to last rib and third to last rib
- Between third to last rib and fourth to last rib
- From second to last rib to fourth to last rib

#### **10th Rib - Depth measurement taken posterior to the trapezius muscle**

- Between 10th and 11th Rib
- Between 11th and 12th Rib
- Between 12th and 13th Rib
- Between 11th and 13th Rib


**Note: Rib numbers are approximate for market weight pigs and start from the same side as the trapezius muscle.**

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# BioSoft Toolbox<sup>®</sup> 4.0

**Biotronics** Inc.  
Ultrasound Solutions

BioSoft Toolbox<sup>®</sup> for Swine 4.0.1.1 © 2007-2020 by Biotronics, Inc.

File Tools Options Help **Exago**

(N) Next Animal

Image Capture Tools  
(S) Stream (F) Freeze  
(C) Capture Image

Scanner Setup  
Depth: 14cm  
Probe Size: 13cm

Image Interpretation Tools  
 (I) IMF  
# of ROI: 2  
ROI Size: 66  
 w/ Depth Assist  
 (R) Area  
 (D) Fat Depth  
(A) Assist With Depths

Sort By Scan Time

Animal List  
16-4  
23-12  
24-9  
24-10  
27-1  
28-1  
28-6  
28-9  
**28-10**  
29-1


14-Oct-20 10:56 AM  
BioSoft Toolbox<sup>®</sup>: 4.0.1.0  
Biotronics DDI for Exago, a licensed product  
Exago Software Version: 01-12

Animal Information

ID	28-10
IMF	0.3
Area	6.83
Fat Depth	0.50
Loin Depth	2.54
Weight	232
Sex	2
Breed	Berkshire
Farrow Date	
Sire	
Comments	

Images Image Data In English Units

Animals: 35 | Animals w/ Images: 35 | Pixel Value: 0 | Hardlock Count: 6144 | File Version: 4.0.1.0 | Session Folder: Z:\BioQProc\Previous Years\2020\AJC Berkshire...

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# BioSoft Toolbox<sup>®</sup> 4.0

**Biotronics** Inc.  
Ultrasound Solutions

BioSoft Toolbox<sup>®</sup> for Swine 4.0.1.1 © 2007-2020 by Biotronics, Inc.

File Tools Options Help

(N) Next Animal

Image Capture Tools  
(S) Stream (F) Freeze  
(C) Capture Image

Scanner Setup  
Depth: 14cm  
Probe Size: 13cm

Image Interpretation Tools  
 (I) IMF  
# of ROI: 2  
ROI Size: 66  
 w/ Depth Assist  
 (R) Area  
 (D) Fat Depth  
(A) Assist With Depths

Sort By Scan Time

Animal List  
16-4  
23-12  
24-9  
24-10  
27-1  
28-1  
28-6  
28-9  
28-10  
29-1

14-Oct-20 10:56 AM  
BioSoft Toolbox<sup>®</sup>: 4.0.1.0  
Biotronics DDI for Exago, a licensed product  
Exago Software Version: 01-12

Animal Information  
ID: 28-10  
IMF: 0.3  
Area: 7.16  
Fat Depth: 0.50  
Loin Depth:  
Weight: 232  
Sex: 2  
Breed: Berkshire  
Farrow Date:  
Sire:  
Comments:

Images Image Data In English Units  
1 2 3 4 5

Animals: 35 | Animals w/ Images: 35 | Pixel Value: 0 | Hardlock Count: 6145 | File Version: 4.0.1.0 | Session Folder: Z:\BioQProc\Previous Years\2020\AJC Berkshire...

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IMV TECHNOLOGIES. FOR LIFE.

Thanks for your attention!!